

stryker[®]

Orthopaedics



Triathlon TS Knee System

Product Catalog

Table of Contents

Femoral Components

Total Stabilizer

Tibial Baseplates

Total Stabilizer

Tibial Inserts

Total Stabilizer

Patellas

Symmetric

- Cemented

Asymmetric

- Cemented & Cementless

Stems

Total Stabilizer

Stem Extenders

Total Stabilizer

Offset Adapters

Total Stabilizer

Augments

Total Stabilizer

Instruments

Total Stabilizer





Orthopaedics

Femoral Components


**Femoral
Components**

stryker[®]

Orthopaedics

Total Stabilizer





Size	A/P Width (mm)	M/L Width (mm)	Resected A/P (mm)	Box Height (mm)
# 1	53	59	37.34	23mm
# 2	56	62	39.95	23mm
# 3	59	65	42.49	23mm
# 4	62	68	45.26	23mm
# 5	65	71	48.44	23mm
# 6	68	74	51.44	23mm
# 7	71	77	53.47	23mm
# 8	75	80	57.73	23mm



Catalog #	Description	Size
5512-F-101	TS Femoral Component	# 1 Left
5512-F-201	TS Femoral Component	# 2 Left
5512-F-301	TS Femoral Component	# 3 Left
5512-F-401	TS Femoral Component	# 4 Left
5512-F-501	TS Femoral Component	# 5 Left
5512-F-601	TS Femoral Component	# 6 Left
5512-F-701	TS Femoral Component	# 7 Left
5512-F-801	TS Femoral Component	# 8 Left

Catalog #	Description	Size
5512-F-102	TS Femoral Component	# 1 Right
5512-F-202	TS Femoral Component	# 2 Right
5512-F-302	TS Femoral Component	# 3 Right
5512-F-402	TS Femoral Component	# 4 Right
5512-F-502	TS Femoral Component	# 5 Right
5512-F-602	TS Femoral Component	# 6 Right
5512-F-702	TS Femoral Component	# 7 Right
5512-F-802	TS Femoral Component	# 8 Right



Orthopaedics

Tibial Baseplates



**Tibial
Baseplates**

stryker[®]

Orthopaedics

Total Stabilizer





Size	A/P Width (mm)	M/L Width (mm)
# 1	40	61
# 2	42	64
# 3	44	67
# 4	46	70
# 5	49	74
# 6	52	77
# 7	56	80
# 8	59	85



Catalog #	Description	Size
5521-B-100	TS Universal Baseplate	# 1
5521-B-200	TS Universal Baseplate	# 2
5521-B-300	TS Universal Baseplate	# 3
5521-B-400	TS Universal Baseplate	# 4
5521-B-500	TS Universal Baseplate	# 5
5521-B-600	TS Universal Baseplate	# 6
5521-B-700	TS Universal Baseplate	# 7
5521-B-800	TS Universal Baseplate	# 8



Orthopaedics

Tibial Inserts


**Tibial
Inserts**

stryker[®]

Orthopaedics

Total Stabilizer





Size (mm)	1-8
Thickness (mm)	9, 11, 13, 16, 19, 22, 25, 28, 31
Post Height (mm)	25.6
*Jump Height (mm)	18
Post Width (mm)	15.6
Varus/Valgus Constraint	± 2
Internal/External Constraint	± 7 at full extension ± 3 at 90°extension
Posterior Slope	0
Cam Engagement	40°

*The distance the femoral component must travel to clear the post in full flexion




Catalog #	Description	Size
5533-G-109	TS Tibial Insert - X3 Poly	# 1 - 9mm
5533-G-111	TS Tibial Insert - X3 Poly	# 1 - 11mm
5533-G-113	TS Tibial Insert - X3 Poly	# 1 - 13mm
5533-G-116	TS Tibial Insert - X3 Poly	# 1 - 16mm
5533-G-119	TS Tibial Insert - X3 Poly	# 1 - 19mm
5533-G-122	TS Tibial Insert - X3 Poly	# 1 - 22mm
5533-G-125	TS Tibial Insert - X3 Poly	# 1 - 25mm
5533-G-128	TS Tibial Insert - X3 Poly	# 1 - 28mm
5533-G-131	TS Tibial Insert - X3 Poly	# 1 - 31mm
5533-G-209	TS Tibial Insert - X3 Poly	# 2 - 9mm
5533-G-211	TS Tibial Insert - X3 Poly	# 2 - 11mm
5533-G-213	TS Tibial Insert - X3 Poly	# 2 - 13mm
5533-G-216	TS Tibial Insert - X3 Poly	# 2 - 16mm
5533-G-219	TS Tibial Insert - X3 Poly	# 2 - 19mm
5533-G-222	TS Tibial Insert - X3 Poly	# 2 - 22mm
5533-G-225	TS Tibial Insert - X3 Poly	# 2 - 25mm
5533-G-228	TS Tibial Insert - X3 Poly	# 2 - 28mm
5533-G-231	TS Tibial Insert - X3 Poly	# 2 - 31mm

Catalog #	Description	Size
5533-G-309	TS Tibial Insert - X3 Poly	# 3 - 9mm
5533-G-311	TS Tibial Insert - X3 Poly	# 3 - 11mm
5533-G-313	TS Tibial Insert - X3 Poly	# 3 - 13mm
5533-G-316	TS Tibial Insert - X3 Poly	# 3 - 16mm
5533-G-319	TS Tibial Insert - X3 Poly	# 3 - 19mm
5533-G-322	TS Tibial Insert - X3 Poly	# 3 - 22mm
5533-G-325	TS Tibial Insert - X3 Poly	# 3 - 25mm
5533-G-328	TS Tibial Insert - X3 Poly	# 3 - 28mm
5533-G-331	TS Tibial Insert - X3 Poly	# 3 - 31mm
5533-G-409	TS Tibial Insert - X3 Poly	# 4 - 9mm
5533-G-411	TS Tibial Insert - X3 Poly	# 4 - 11mm
5533-G-413	TS Tibial Insert - X3 Poly	# 4 - 13mm
5533-G-416	TS Tibial Insert - X3 Poly	# 4 - 16mm
5533-G-419	TS Tibial Insert - X3 Poly	# 4 - 19mm
5533-G-422	TS Tibial Insert - X3 Poly	# 4 - 22mm
5533-G-425	TS Tibial Insert - X3 Poly	# 4 - 25mm
5533-G-428	TS Tibial Insert - X3 Poly	# 4 - 28mm
5533-G-431	TS Tibial Insert - X3 Poly	# 4 - 31mm



Catalog #	Description	Size
5533-G-509	TS Tibial Insert - X3 Poly	# 5 - 9mm
5533-G-511	TS Tibial Insert - X3 Poly	# 5 - 11mm
5533-G-513	TS Tibial Insert - X3 Poly	# 5 - 13mm
5533-G-516	TS Tibial Insert - X3 Poly	# 5 - 16mm
5533-G-519	TS Tibial Insert - X3 Poly	# 5 - 19mm
5533-G-522	TS Tibial Insert - X3 Poly	# 5 - 22mm
5533-G-525	TS Tibial Insert - X3 Poly	# 5 - 25mm
5533-G-528	TS Tibial Insert - X3 Poly	# 5 - 28mm
5533-G-531	TS Tibial Insert - X3 Poly	# 5 - 31mm
5533-G-609	TS Tibial Insert - X3 Poly	# 6 - 9mm
5533-G-611	TS Tibial Insert - X3 Poly	# 6 - 11mm
5533-G-613	TS Tibial Insert - X3 Poly	# 6 - 13mm
5533-G-616	TS Tibial Insert - X3 Poly	# 6 - 16mm
5533-G-619	TS Tibial Insert - X3 Poly	# 6 - 19mm
5533-G-622	TS Tibial Insert - X3 Poly	# 6 - 22mm
5533-G-625	TS Tibial Insert - X3 Poly	# 6 - 25mm
5533-G-628	TS Tibial Insert - X3 Poly	# 6 - 28mm
5533-G-631	TS Tibial Insert - X3 Poly	# 6 - 31mm

Catalog #	Description	Size
5533-G-709	TS Tibial Insert - X3 Poly	# 7 - 9mm
5533-G-711	TS Tibial Insert - X3 Poly	# 7 - 11mm
5533-G-713	TS Tibial Insert - X3 Poly	# 7 - 13mm
5533-G-716	TS Tibial Insert - X3 Poly	# 7 - 16mm
5533-G-719	TS Tibial Insert - X3 Poly	# 7 - 19mm
5533-G-722	TS Tibial Insert - X3 Poly	# 7 - 22mm
5533-G-725	TS Tibial Insert - X3 Poly	# 7 - 25mm
5533-G-728	TS Tibial Insert - X3 Poly	# 7 - 28mm
5533-G-731	TS Tibial Insert - X3 Poly	# 7 - 31mm
5533-G-809	TS Tibial Insert - X3 Poly	# 8 - 9mm
5533-G-811	TS Tibial Insert - X3 Poly	# 8 - 11mm
5533-G-813	TS Tibial Insert - X3 Poly	# 8 - 13mm
5533-G-816	TS Tibial Insert - X3 Poly	# 8 - 16mm
5533-G-819	TS Tibial Insert - X3 Poly	# 8 - 19mm
5533-G-822	TS Tibial Insert - X3 Poly	# 8 - 22mm
5533-G-825	TS Tibial Insert - X3 Poly	# 8 - 25mm
5533-G-828	TS Tibial Insert - X3 Poly	# 8 - 28mm
5533-G-831	TS Tibial Insert - X3 Poly	# 8 - 31mm




Size (mm)	Diameter (mm)	Thickness (mm)	Peg Length (mm)	Peg Diameter (mm)
S27x8	27	8	5	5-3/4
S29x8	29	8	5	5-3/4
S31x9	31	9	5	5-3/4
S33x9	33	9	5	5-3/4
S36x10	36	10	5	5-3/4
S39x11	39	11	5	5-3/4




Catalog #	Description	Size
5550-L-278	Triathlon Symmetric All Poly Patella	S27x8
5550-L-298	Triathlon Symmetric All Poly Patella	S29x8
5550-L-319	Triathlon Symmetric All Poly Patella	S31x9
5550-L-339	Triathlon Symmetric All Poly Patella	S33x9
5550-L-360	Triathlon Symmetric All Poly Patella	S36x10
5550-L-391	Triathlon Symmetric All Poly Patella	S39x11

Triathlon Asymmetric
All Poly Patella



Triathlon Asymmetric
Press Fit Patella (Metal Backed)



Size (mm)	S/I Width (mm)	M/L Width (mm)	Thickness (mm)	Peg Length (mm)	Peg Diameter (mm)	Peg Diameter (mm)
A29x9	29	33	9	5	5 3/4	N/A
A32x10	32	36	10	5	5 3/4	6
A35x10	35	39	10	5	5 3/4	6
A38x11	38	42	11	5	5 3/4	6
A40x11	40	44	11	5	5 3/4	6



Catalog #	Description	Size
-----------	-------------	------

Triathlon Asymmetric All Poly Patella

5551-L-299	Triathlon Asymmetric All Poly Patella	S27x8
5551-L-320	Triathlon Asymmetric All Poly Patella	S29x8
5551-L-350	Triathlon Asymmetric All Poly Patella	S31x9
5551-L-381	Triathlon Asymmetric All Poly Patella	S33x9
5551-L-401	Triathlon Asymmetric All Poly Patella	S36x10

Triathlon Asymmetric Metal-Backed Patella

5554-L-320	Triathlon Asymmetric Metal-Backed Patella	A32x10
5554-L-350	Triathlon Asymmetric Metal-Backed Patella	A35x10
5554-L-381	Triathlon Asymmetric Metal-Backed Patella	A38x11
5554-L-401	Triathlon Asymmetric Metal-Backed Patella	A40x11

stryker[®]

Orthopaedics

Stems


Stems

stryker[®]

Orthopaedics

Total Stabilizer






Stem Type	Diameter (mm)	Length (mm)
Cemented	9, 12, 15	50, 100, 150
Cementless	10-25 (1mm increments)	100, 150
Stem Extenders	16	25, 50



Catalog #	Description	Size
5560-S-109	Triathlon Cemented Stem and End Cap	9mm x 50mm
5560-S-209	Triathlon Cemented Stem and End Cap	9mm x 100mm
5560-S-309	Triathlon Cemented Stem and End Cap	9mm x 150mm
5560-S-112	Triathlon Cemented Stem and End Cap	12mm x 50mm
5560-S-212	Triathlon Cemented Stem and End Cap	12mm x 100mm
5560-S-312	Triathlon Cemented Stem and End Cap	12mm x 150mm
5560-S-115	Triathlon Cemented Stem and End Cap	15mm x 50mm
5560-S-215	Triathlon Cemented Stem and End Cap	15mm x 100mm
5560-S-315	Triathlon Cemented Stem and End Cap	15mm x 150mm



Stem Type	Diameter (mm)	Length (mm)
Cemented	9, 12, 15	50, 100, 150
Cementless	10-25	100, 150
Stem Extenders	16	25, 50



Catalog #	Description	Size
5565-S-010	Triathlon Cementless Stem & End Cap	10mm x 100mm
5565-S-011	Triathlon Cementless Stem & End Cap	11mm x 100mm
5565-S-012	Triathlon Cementless Stem & End Cap	12mm x 100mm
5565-S-013	Triathlon Cementless Stem & End Cap	13mm x 100mm
5565-S-014	Triathlon Cementless Stem & End Cap	14mm x 100mm
5565-S-015	Triathlon Cementless Stem & End Cap	15mm x 100mm
5565-S-016	Triathlon Cementless Stem & End Cap	16mm x 100mm
5565-S-017	Triathlon Cementless Stem & End Cap	17mm x 100mm
5565-S-018	Triathlon Cementless Stem & End Cap	18mm x 100mm
5565-S-019	Triathlon Cementless Stem & End Cap	19mm x 100mm
5565-S-020	Triathlon Cementless Stem & End Cap	20mm x 100mm
5565-S-021	Triathlon Cementless Stem & End Cap	21mm x 100mm
5565-S-022	Triathlon Cementless Stem & End Cap	22mm x 100mm
5565-S-023	Triathlon Cementless Stem & End Cap	23mm x 100mm
5565-S-024	Triathlon Cementless Stem & End Cap	24mm x 100mm
5565-S-025	Triathlon Cementless Stem & End Cap	25mm x 100mm

Catalog #	Description	Size
5566-S-010	Triathlon Cementless Stem & End Cap	10mm x 150mm
5566-S-011	Triathlon Cementless Stem & End Cap	11mm x 150mm
5566-S-012	Triathlon Cementless Stem & End Cap	12mm x 150mm
5566-S-013	Triathlon Cementless Stem & End Cap	13mm x 150mm
5566-S-014	Triathlon Cementless Stem & End Cap	14mm x 150mm
5566-S-015	Triathlon Cementless Stem & End Cap	15mm x 150mm
5566-S-016	Triathlon Cementless Stem & End Cap	16mm x 150mm
5566-S-017	Triathlon Cementless Stem & End Cap	17mm x 150mm
5566-S-018	Triathlon Cementless Stem & End Cap	18mm x 150mm
5566-S-019	Triathlon Cementless Stem & End Cap	19mm x 150mm
5566-S-020	Triathlon Cementless Stem & End Cap	20mm x 150mm
5566-S-021	Triathlon Cementless Stem & End Cap	21mm x 150mm
5566-S-022	Triathlon Cementless Stem & End Cap	22mm x 150mm
5566-S-023	Triathlon Cementless Stem & End Cap	23mm x 150mm
5566-S-024	Triathlon Cementless Stem & End Cap	24mm x 150mm
5566-S-025	Triathlon Cementless Stem & End Cap	25mm x 150mm



Orthopaedics

Stem Extenders/ Offset Adapters


**Stem Extenders/
Offset Adapters**

stryker[®]

Orthopaedics

Total Stabilizer






Stem Type	Diameter (mm)	Length (mm)
Cemented	9, 12, 15	50, 100, 150
Cementless	10-25	100, 150
Stem Extenders	16	25, 50



Catalog #	Description	Size
5571-S-025	Triathlon Stem Extender	25mm
5571-S-050	Triathlon Stem Extender	50mm



Size (mm)	Stem Length (mm)
2, 4, 6, 8	25



Catalog #	Description	Size
5570-S-020	Triathlon TS Offset Adapter	2mm
5570-S-040	Triathlon TS Offset Adapter	4mm
5570-S-060	Triathlon TS Offset Adapter	6mm
5570-S-080	Triathlon TS Offset Adapter	8mm

stryker[®]

Orthopaedics

Augments


Augments

Total Stabilizer

stryker®

Orthopaedics






LEFT 5mm #4




RIGHT 5mm #4




LEFT 10mm #4




RIGHT 10mm #4



LEFT 15mm #4



RIGHT 15mm #4



Size (mm)	Left Distal (mm)	Right Distal (mm)	Half Posterior (mm)
1-8	5, 10, 15	5, 10, 15	5, 10





Catalog #	Description	Size
5540-A-101	Triathlon Femoral Distal Augment - 5mm	# 1 Left
5540-A-201	Triathlon Femoral Distal Augment - 5mm	# 2 Left
5540-A-301	Triathlon Femoral Distal Augment - 5mm	# 3 Left
5540-A-401	Triathlon Femoral Distal Augment - 5mm	# 4 Left
5540-A-501	Triathlon Femoral Distal Augment - 5mm	# 5 Left
5540-A-601	Triathlon Femoral Distal Augment - 5mm	# 6 Left
5540-A-701	Triathlon Femoral Distal Augment - 5mm	# 7 Left
5540-A-801	Triathlon Femoral Distal Augment - 5mm	# 8 Left
5540-A-102	Triathlon Femoral Distal Augment - 5mm	# 1 Right
5540-A-202	Triathlon Femoral Distal Augment - 5mm	# 2 Right
5540-A-302	Triathlon Femoral Distal Augment - 5mm	# 3 Right
5540-A-402	Triathlon Femoral Distal Augment - 5mm	# 4 Right
5540-A-502	Triathlon Femoral Distal Augment - 5mm	# 5 Right
5540-A-602	Triathlon Femoral Distal Augment - 5mm	# 6 Right
5540-A-702	Triathlon Femoral Distal Augment - 5mm	# 7 Right
5540-A-802	Triathlon Femoral Distal Augment - 5mm	# 8 Right

Catalog #	Description	Size
5541-A-101	Triathlon Femoral Distal Augment - 10mm	# 1 Left
5541-A-201	Triathlon Femoral Distal Augment - 10mm	# 2 Left
5541-A-301	Triathlon Femoral Distal Augment - 10mm	# 3 Left
5541-A-401	Triathlon Femoral Distal Augment - 10mm	# 4 Left
5541-A-501	Triathlon Femoral Distal Augment - 10mm	# 5 Left
5541-A-601	Triathlon Femoral Distal Augment - 10mm	# 6 Left
5541-A-701	Triathlon Femoral Distal Augment - 10mm	# 7 Left
5541-A-801	Triathlon Femoral Distal Augment - 10mm	# 8 Left
5541-A-102	Triathlon Femoral Distal Augment - 10mm	# 1 Right
5541-A-202	Triathlon Femoral Distal Augment - 10mm	# 2 Right
5541-A-302	Triathlon Femoral Distal Augment - 10mm	# 3 Right
5541-A-402	Triathlon Femoral Distal Augment - 10mm	# 4 Right
5541-A-502	Triathlon Femoral Distal Augment - 10mm	# 5 Right
5541-A-602	Triathlon Femoral Distal Augment - 10mm	# 6 Right
5541-A-702	Triathlon Femoral Distal Augment - 10mm	# 7 Right
5541-A-802	Triathlon Femoral Distal Augment - 10mm	# 8 Right





Catalog #	Description	Size
5542-A-101	Triathlon Femoral Distal Augment - 15mm	# 1 Left
5542-A-201	Triathlon Femoral Distal Augment - 15mm	# 2 Left
5542-A-301	Triathlon Femoral Distal Augment - 15mm	# 3 Left
5542-A-401	Triathlon Femoral Distal Augment - 15mm	# 4 Left
5542-A-501	Triathlon Femoral Distal Augment - 15mm	# 5 Left
5542-A-601	Triathlon Femoral Distal Augment - 15mm	# 6 Left
5542-A-701	Triathlon Femoral Distal Augment - 15mm	# 7 Left
5542-A-801	Triathlon Femoral Distal Augment - 15mm	# 8 Left
5542-A-102	Triathlon Femoral Distal Augment - 15mm	# 1 Right
5542-A-202	Triathlon Femoral Distal Augment - 15mm	# 2 Right
5542-A-302	Triathlon Femoral Distal Augment - 15mm	# 3 Right
5542-A-402	Triathlon Femoral Distal Augment - 15mm	# 4 Right
5542-A-502	Triathlon Femoral Distal Augment - 15mm	# 5 Right
5542-A-602	Triathlon Femoral Distal Augment - 15mm	# 6 Right
5542-A-702	Triathlon Femoral Distal Augment - 15mm	# 7 Right
5542-A-802	Triathlon Femoral Distal Augment - 15mm	# 8 Right

Catalog #	Description	Size
Triathlon Femoral Augment Locking Screw		
5540-A-000	Triathlon Femoral Augment Locking Screw	
5545-A-000	Triathlon Femoral Augment Locking Screw	



LEFT MEDIAL / RIGHT LATERAL



RIGHT MEDIAL / LEFT LATERAL

Size (mm)	Half Block (mm)
1-8	5, 10



Catalog #	Description	Size
5545-A-101	Triathlon Tibial Augment - 5mm	# 1 Left
5545-A-201	Triathlon Tibial Augment - 5mm	# 2 Left
5545-A-301	Triathlon Tibial Augment - 5mm	# 3 Left
5545-A-401	Triathlon Tibial Augment - 5mm	# 4 Left
5545-A-501	Triathlon Tibial Augment - 5mm	# 5 Left
5545-A-601	Triathlon Tibial Augment - 5mm	# 6 Left
5545-A-701	Triathlon Tibial Augment - 5mm	# 7 Left
5545-A-801	Triathlon Tibial Augment - 5mm	# 8 Left
5545-A-102	Triathlon Tibial Augment - 5mm	# 1 Right
5545-A-202	Triathlon Tibial Augment - 5mm	# 2 Right
5545-A-302	Triathlon Tibial Augment - 5mm	# 3 Right
5545-A-402	Triathlon Tibial Augment - 5mm	# 4 Right
5545-A-502	Triathlon Tibial Augment - 5mm	# 5 Right
5545-A-602	Triathlon Tibial Augment - 5mm	# 6 Right
5545-A-702	Triathlon Tibial Augment - 5mm	# 7 Right
5545-A-802	Triathlon Tibial Augment - 5mm	# 8 Right

Catalog #	Description	Size
5546-A-101	Triathlon Tibial Augment - 10mm	# 1 Left
5546-A-201	Triathlon Tibial Augment - 10mm	# 2 Left
5546-A-301	Triathlon Tibial Augment - 10mm	# 3 Left
5546-A-401	Triathlon Tibial Augment - 10mm	# 4 Left
5546-A-501	Triathlon Tibial Augment - 10mm	# 5 Left
5546-A-601	Triathlon Tibial Augment - 10mm	# 6 Left
5546-A-701	Triathlon Tibial Augment - 10mm	# 7 Left
5546-A-801	Triathlon Tibial Augment - 10mm	# 8 Left
5546-A-102	Triathlon Tibial Augment - 10mm	# 1 Right
5546-A-202	Triathlon Tibial Augment - 10mm	# 2 Right
5546-A-302	Triathlon Tibial Augment - 10mm	# 3 Right
5546-A-402	Triathlon Tibial Augment - 10mm	# 4 Right
5546-A-502	Triathlon Tibial Augment - 10mm	# 5 Right
5546-A-602	Triathlon Tibial Augment - 10mm	# 6 Right
5546-A-702	Triathlon Tibial Augment - 10mm	# 7 Right
5546-A-802	Triathlon Tibial Augment - 10mm	# 8 Right

stryker[®]

Orthopaedics

Instruments


Instruments

stryker[®]


Orthopaedics

Total Stabilizer







Catalog #	Description	Qty	Ref. Letter
6541-4-810	Impaction Handle	2	A
6541-4-825	Slip Torque Handle	1	B
6543-7-601	Resection Guide Tower	1	C
6541-4-807	Femoral Impactor Extractor	1	D
6541-4-811	Femoral Impactor	1	E
6541-4-812	Tibial Baseplate Impactor	1	F
6541-4-813	Tibial Insert Impactor	1	G
6541-4-805	Tibial Baseplate Impactor Extractor	1	H
6543-4-516	Stem Extender Shaft	1	I
6541-4-806	Universal Alignment Handle	1	J
6543-4-802	1/8" Universal Hex Driver	1	K
6541-4-516	5/16" IM Rod	1	L
6541-4-602	Universal Alignment Rod	1	M
6543-7-600	Support Arm Assembly	1	N



Catalog #	Description	Qty	Ref. Letter
6543-4-803	Offset Counter Wrench	1	A
6543-4-801	Universal Counter Wrench	1	B
6543-4-818	Universal Torque Wrench	1	C
6541-4-400	Bladerunner	1	D
6543-7-602	Stop Plate	1	E
6543-4-605	Adjustable Spacer Block Augment-5mm	4	F
6543-4-610	Adjustable Spacer Block Augment-10mm	4	G
6543-4-610	Adjustable Spacer Block Augment-15mm	2	H
6541-4-610	Adjustable Spacer Block	1	I
6541-4-804	Headless Pin Extractor	1	J
6541-4-300	Headed Nail Impactor Extractor	1	K
6541-4-803	Slap Hammer	1	L
6541-4-515	Headed Nails - 1 1/2"	2	M
6541-4-575	Headed Nails - 3/4"	2	N



Catalog #	Description	Qty	Ref. Letter
6543-4-527	Boss/Offset Reamer	1	A
6543-7-508	8mm Starter Awl	1	B
6541-4-518	1/8" Drill	1	C
6541-4-800	T-Handle Driver	1	D
6541-4-801	Universal Driver	1	E
6541-4-809	Headless Pin Driver	1	F
6541-4-003	Headless Pins - 3"	4	G
6543-7-509	IM Reamer - 9mm	1	H
6543-7-510	IM Reamer - 10mm	1	I
6543-7-511	IM Reamer - 11mm	1	J
6543-7-512	IM Reamer - 12mm	1	K
6543-7-513	IM Reamer - 13mm	1	L
6543-7-514	IM Reamer - 14mm	1	M
6543-7-515	IM Reamer - 15mm	1	N
5560-T-112	12 x 50mm Cemented Stem Trial	2	O




Catalog #	Description	Qty	Ref. Letter
6543-7-516	IM Reamer - 16mm	1	A
6543-7-517	IM Reamer - 17mm	1	B
6543-7-518	IM Reamer - 18mm	1	C
6543-7-519	IM Reamer - 19mm	1	D
6543-7-520	IM Reamer - 20mm	1	E
6543-7-521	IM Reamer - 21mm	1	F
5560-T-115	15mm x 50mm	2	G
5571-T-025	Triathlon Stem Extender Trial - 25mm	2	H
5571-T-050	Triathlon Stem Extender Trial - 50mm	2	I




Catalog #	Description	Qty	Ref. Letter
5565-T-009	Triathlon Stem Trial, 9 x 100mm	1	A
5565-T-010	Triathlon Stem Trial, 10 x 100mm	1	B
5565-T-011	Triathlon Stem Trial, 11 x 100mm	1	C
5565-T-012	Triathlon Stem Trial, 12 x 100mm	2	D
5565-T-013	Triathlon Stem Trial, 13 x 100mm	2	E
5565-T-014	Triathlon Stem Trial, 14 x 100mm	2	F
5565-T-015	Triathlon Stem Trial, 15 x 100mm	2	G
5565-T-016	Triathlon Stem Trial, 16 x 100mm	2	H
5565-T-017	Triathlon Stem Trial, 17 x 100mm	2	I
5565-T-018	Triathlon Stem Trial, 18 x 100mm	2	J
5565-T-019	Triathlon Stem Trial, 19 x 100mm	2	K
5565-T-020	Triathlon Stem Trial, 20 x 100mm	2	L
5565-T-021	Triathlon Stem Trial, 21 x 100mm	2	M

Catalog #	Description	Qty	Ref. Letter
5566-T-009	Triathlon Stem Trial, 9 x 150mm	1	N
5566-T-010	Triathlon Stem Trial, 10 x 150mm	1	O
5566-T-011	Triathlon Stem Trial, 11 x 150mm	1	P
5566-T-012	Triathlon Stem Trial, 12 x 150mm	2	Q
5566-T-013	Triathlon Stem Trial, 13 x 150mm	2	R
5566-T-014	Triathlon Stem Trial, 14 x 150mm	2	S
5566-T-015	Triathlon Stem Trial, 15 x 150mm	2	T
5566-T-016	Triathlon Stem Trial, 16 x 150mm	2	U
5566-T-017	Triathlon Stem Trial, 17 x 150mm	2	V
5566-T-018	Triathlon Stem Trial, 18 x 150mm	2	W
5566-T-019	Triathlon Stem Trial, 19 x 150mm	2	X
5566-T-020	Triathlon Stem Trial, 20 x 150mm	2	Y
5566-T-021	Triathlon Stem Trial, 21 x 150mm	2	Z




Catalog #	Description	Qty	Ref. Letter
6543-1-005	Distal Spacer - 5mm	1	A
6543-1-010	Distal Spacer - 10mm	1	B
6543-1-015	Distal Spacer - 15mm	1	C
6543-1-600	Femoral Offset Bushing	1	D
6543-1-603	Size 1-8 Femoral Sizing Templates	1	E
6543-1-703	#3 All-in-One Cutting Block	1	F
6543-1-704	#4 All-in-One Cutting Block	1	G
6543-1-705	#5 All-in-One Cutting Block	1	H
6543-1-706	#6 All-in-One Cutting Block	1	I
6543-1-710	Revision Box Cutting Guide	1	J
6543-1-721	Revision Distal Resection Guide	1	K
6543-1-750	Femoral Boss Preparation Guide	1	L
6543-1-751	Femoral Boss Bushing	1	M
6543-4-400	Joint Line Ruler	1	N




Catalog #	Description	Qty	Ref. Letter
5512-T-301	TS Femoral Trials	1	A
5512-T-302	TS Femoral Trials	1	B
5512-T-401	TS Femoral Trials	1	C
5512-T-402	TS Femoral Trials	1	D
5512-T-501	TS Femoral Trials	1	G
5512-T-501	TS Femoral Trials	1	H
5512-T-601	TS Femoral Trials	1	E
5512-T-602	TS Femoral Trials	1	F
5570-T-020	Triathlon Offset Adaptor Trial - 2mm	1	I
5570-T-040	Triathlon Offset Adaptor Trial - 4mm	1	J
5570-T-060	Triathlon Offset Adaptor Trial - 6mm	1	K
5570-T-080	Triathlon Offset Adaptor Trial - 8mm	1	L
5540-T-300	Triathlon Femoral Distal Augment Trial, 5mm - # 3	2	M
5540-T-400	Triathlon Femoral Distal Augment Trial, 5mm - # 4	2	N
5540-T-500	Triathlon Femoral Distal Augment Trial, 5mm - # 5	2	P
5540-T-600	Triathlon Femoral Distal Augment Trial, 5mm - # 6	2	O
5541-T-300	Triathlon Femoral Distal Augment Trial, 10mm - # 3	2	Q

Catalog #	Description	Qty	Ref. Letter
5541-T-400	Triathlon Femoral Distal Augment Trial, 10mm - # 4	2	R
5541-T-500	Triathlon Femoral Distal Augment Trial, 10mm - # 5	2	T
5541-T-600	Triathlon Femoral Distal Augment Trial, 10mm - # 6	2	S
5542-T-300	Triathlon Femoral Distal Augment Trial, 15mm - # 3	2	U
5542-T-400	Triathlon Femoral Distal Augment Trial, 15mm - # 4	2	V
5542-T-500	Triathlon Femoral Distal Augment Trial, 15mm - # 5	2	X
5542-T-600	Triathlon Femoral Distal Augment Trial, 15mm - # 6	2	W
5543-T-300	Triathlon Femoral Posterior Augment Trial, 5mm - # 3	2	Y
5543-T-400	Triathlon Femoral Posterior Augment Trial, 5mm - # 4	2	Z
5543-T-500	Triathlon Femoral Posterior Augment Trial, 5mm - # 5	2	AA
5543-T-600	Triathlon Femoral Posterior Augment Trial, 5mm - # 6	2	AB
5544-T-300	Triathlon Femoral Posterior Augment Trial, 10mm - # 3	2	AC
5544-T-400	Triathlon Femoral Posterior Augment Trial, 10mm - # 4	2	AD
5544-T-500	Triathlon Femoral Posterior Augment Trial, 10mm - # 5	2	AE
5544-T-600	Triathlon Femoral Posterior Augment Trial, 10mm - # 6	2	AF




Catalog #	Description	Qty	Ref. Letter
6543-2-600	Tibial Offset Bushing	1	A
6543-2-601	Tibial Offset Bushing Guide	1	B
6543-6-700	Revision Tibial Resection Guide L - Slotted	1	C
6543-6-701	Revision Tibial Resection Guide R - Slotted	1	D
6541-2-013	Size 1-3 Keel Punch	1	E
6541-2-046	Size 4-6 Keel Punch	1	F
6541-2-603	#3 Universal Tibial Template	1	G
6541-2-604	#4 Universal Tibial Template	1	H
6541-2-605	#5 Universal Tibial Template	1	I
6541-2-606	#6 Universal Tibial Template	1	J
6541-2-713	Size 1-3 Keel Punch Guide	1	K
6541-2-748	Size 4-8 Keel Punch Guide	1	L
6541-2-807	Tibial Alignment Handle	1	M
6543-2-703	Tibial Resection Guide Link	1	N




Catalog #	Description	Qty	Ref. Letter
5521-T-300	TS Baseplate Trials # 3 - 9mm	1	A
5521-T-400	TS Baseplate Trials # 4 - 9mm	1	B
5521-T-500	TS Baseplate Trials # 5 - 9mm	1	C
5521-T-600	TS Baseplate Trials # 6 - 9mm	1	D
5545-T-301	Tibial Augment Trial #3 LM/RL - 5mm	1	E
5545-T-302	Tibial Augment Trial #3 RM/LL - 5mm	1	H
5545-T-401	Tibial Augment Trial #4 LM/RL - 5mm	1	I
5545-T-402	Tibial Augment Trial #4 RM/LL - 5mm	1	L
5545-T-501	Tibial Augment Trial #5 LM/RL - 5mm	1	M
5545-T-502	Tibial Augment Trial #5 RM/LL - 5mm	1	P
5545-T-601	Tibial Augment Trial #6 LM/RL - 5mm	1	Q
5545-T-602	Tibial Augment Trial #6 RM/LL - 5mm	1	T

Catalog #	Description	Qty	Ref. Letter
5546-T-301	Tibial Augment Trial #3 LM/RL - 10mm	1	F
5546-T-302	Tibial Augment Trial #3 RM/LL - 10mm	1	G
5546-T-401	Tibial Augment Trial #4 LM/RL - 10mm	1	J
5546-T-402	Tibial Augment Trial #4 RM/LL - 10mm	1	K
5546-T-501	Tibial Augment Trial #5 LM/RL - 10mm	1	N
5546-T-502	Tibial Augment Trial #5 RM/LL - 10mm	1	O
5546-T-601	Tibial Augment Trial #6 LM/RL - 10mm	1	R
5546-T-602	Tibial Augment Trial #6 RM/LL - 10mm	1	S
5570-T-020	Triathlon Offset Adaptor Trial - 2mm	1	U
5570-T-040	Triathlon Offset Adaptor Trial - 4mm	1	V
5570-T-060	Triathlon Offset Adaptor Trial - 6mm	1	W
5570-T-080	Triathlon Offset Adaptor Trial - 8mm	1	X




Catalog #	Description	Qty	Ref. Letter
5537-T-309	Triathlon TS+ Tibial Insert Trial #3 - 9mm	1	A
5537-T-311	Triathlon TS+ Tibial Insert Trial #3 - 11mm	1	B
5537-T-313	Triathlon TS+ Tibial Insert Trial #3 - 13mm	1	C
5537-T-316	Triathlon TS+ Tibial Insert Trial #3 - 16mm	1	D
5537-T-319	Triathlon TS+ Tibial Insert Trial #3 - 19mm	1	E
5537-T-322	Triathlon TS+ Tibial Insert Trial #3 - 22mm	1	F
5537-T-325	Triathlon TS+ Tibial Insert Trial #3 - 25mm	1	G
5537-T-328	Triathlon TS+ Tibial Insert Trial #3 - 28mm	1	H
5537-T-331	Triathlon TS+ Tibial Insert Trial #3 - 31mm	1	I
5537-T-409	Triathlon TS+ Tibial Insert Trial #4 - 9mm	1	J
5537-T-411	Triathlon TS+ Tibial Insert Trial #4 - 11mm	1	K
5537-T-413	Triathlon TS+ Tibial Insert Trial #4 - 13mm	1	L
5537-T-416	Triathlon TS+ Tibial Insert Trial #4 - 16mm	1	M
5537-T-419	Triathlon TS+ Tibial Insert Trial #4 - 19mm	1	N
5537-T-422	Triathlon TS+ Tibial Insert Trial #4 - 22mm	1	O
5537-T-425	Triathlon TS+ Tibial Insert Trial #4 - 25mm	1	P
5537-T-428	Triathlon TS+ Tibial Insert Trial #4 - 28mm	1	Q
5537-T-431	Triathlon TS+ Tibial Insert Trial #4 - 31mm	1	R

Catalog #	Description	Qty	Ref. Letter
5537-T-509	Triathlon TS+ Tibial Insert Trial #5 - 9mm	1	S
5537-T-511	Triathlon TS+ Tibial Insert Trial #5 - 11mm	1	T
5537-T-513	Triathlon TS+ Tibial Insert Trial #5 - 13mm	1	U
5537-T-516	Triathlon TS+ Tibial Insert Trial #5 - 16mm	1	V
5537-T-519	Triathlon TS+ Tibial Insert Trial #5 - 19mm	1	W
5537-T-522	Triathlon TS+ Tibial Insert Trial #5 - 22mm	1	X
5537-T-525	Triathlon TS+ Tibial Insert Trial #5 - 25mm	1	Y
5537-T-528	Triathlon TS+ Tibial Insert Trial #5 - 28mm	1	Z
5537-T-531	Triathlon TS+ Tibial Insert Trial #5 - 31mm	1	AA
5537-T-609	Triathlon TS+ Tibial Insert Trial #6 - 9mm	1	AB
5537-T-611	Triathlon TS+ Tibial Insert Trial #6 - 11mm	1	AC
5537-T-613	Triathlon TS+ Tibial Insert Trial #6 - 13mm	1	AD
5537-T-616	Triathlon TS+ Tibial Insert Trial #6 - 16mm	1	AE
5537-T-619	Triathlon TS+ Tibial Insert Trial #6 - 19mm	1	AF
5537-T-622	Triathlon TS+ Tibial Insert Trial #6 - 22mm	1	AG
5537-T-625	Triathlon TS+ Tibial Insert Trial #6 - 25mm	1	AH
5537-T-628	Triathlon TS+ Tibial Insert Trial #6 - 28mm	1	AI
5537-T-631	Triathlon TS+ Tibial Insert Trial #6 - 31mm	1	AJ




Catalog #	Description	Qty	Ref. Letter
6543-1-702	#2 All-in-One Cutting Block	1	A
6543-1-707	#7 All-in-One Cutting Block	1	B
6541-2-078	Size 7-8 Keel Punch	1	C
6541-2-602	#2 Universal Tibial Template	1	D
6541-2-607	#7 Universal Tibial Template	1	E
5521-T-200	#2 Baseplate Trial	1	F
5521-T-700	#7 Baseplate Trial	1	G
5512-T-201	#2 Femoral Trial Left	1	H
5512-T-202	#2 Femoral Trial Right	1	I
5512-T-801	#7 Femoral Trial Left	1	J
5512-T-802	#7 Femoral Trial Right	1	K
5540-T-200	Triathlon Femoral Distal Augment Trial, 5mm - # 2	2	L
5540-T-700	Triathlon Femoral Distal Augment Trial, 5mm - # 7	2	M
5541-T-200	Triathlon Femoral Distal Augment Trial, 10mm - # 2	2	N
5541-T-700	Triathlon Femoral Distal Augment Trial, 10mm - # 7	2	O

Catalog #	Description	Qty	Ref. Letter
5542-T-200	Triathlon Femoral Distal Augment Trial, 15mm - # 2	2	P
5542-T-700	Triathlon Femoral Distal Augment Trial, 15mm - # 7	2	Q
5543-T-200	Triathlon Femoral Posterior Augment Trial, 5mm - # 2	2	R
5543-T-700	Triathlon Femoral Posterior Augment Trial, 5mm - # 7	2	S
5544-T-200	Triathlon Femoral Posterior Augment Trial, 10mm - # 2	2	T
5544-T-700	Triathlon Femoral Posterior Augment Trial, 10mm - # 7	2	U
5545-T-701	Tibial Augment Trial #7 LM/RL - 5mm	1	V
5545-T-702	Tibial Augment Trial #7 RM/LL - 5mm	1	W
5546-T-701	Tibial Augment Trial #7 LM/RL - 10mm	1	X
5546-T-702	Tibial Augment Trial #7 RM/LL - 10mm	1	Y
5545-T-201	Tibial Augment Trial #2 LM/RL - 5mm	1	Z
5545-T-202	Tibial Augment Trial #2 RM/LL - 5mm	1	AA
5546-T-201	Tibial Augment Trial #2 LM/RL - 10mm	1	AB
5546-T-202	Tibial Augment Trial #2 RM/LL - 10mm	1	AC




Catalog #	Description	Qty	Ref. Letter
6543-1-701	#1 All-in-One Cutting Block	1	A
6543-1-708	#8 All-in-One Cutting Block	1	B
6541-2-601	#1 Universal Tibial Template	1	C
6541-2-608	#8 Universal Tibial Template	1	D
5521-T-800	#8 Baseplate Trial	1	E
5521-T-100	#1 Baseplate Trial	1	F
5512-T-101	#1 Femoral Trial Left	1	G
5512-T-102	#1 Femoral Trial Right	1	H
5512-T-801	#8 Femoral Trial Left	1	I
5512-T-802	#8 Femoral Trial Right	1	J
5540-T-100	Triathlon Femoral Distal Augment Trial, 5mm - #1	2	K
5540-T-800	Triathlon Femoral Distal Augment Trial, 5mm - #8	2	L
5541-T-100	Triathlon Femoral Distal Augment Trial, 10mm - #1	2	M
5541-T-800	Triathlon Femoral Distal Augment Trial, 10mm - #8	2	N

Catalog #	Description	Qty	Ref. Letter
5542-T-100	Triathlon Femoral Distal Augment Trial, 15mm - #1	2	O
5542-T-800	Triathlon Femoral Distal Augment Trial, 15mm - #8	2	P
5543-T-100	Triathlon Femoral Posterior Augment Trial, 5mm - #1	2	Q
5543-T-800	Triathlon Femoral Posterior Augment Trial, 5mm - #8	2	R
5544-T-100	Triathlon Femoral Posterior Augment Trial, 10mm - #1	2	S
5544-T-800	Triathlon Femoral Posterior Augment Trial, 10mm - #8	2	T
5545-T-101	Tibial Augment Trial #1 LM/RL - 5mm	1	U
5545-T-102	Tibial Augment Trial #1 RM/LL - 5mm	1	V
5546-T-101	Tibial Augment Trial #1 LM/RL - 10mm	1	W
5546-T-102	Tibial Augment Trial #1 RM/LL - 10mm	1	X
5545-T-801	Tibial Augment Trial #8 LM/RL - 5mm	1	Y
5545-T-802	Tibial Augment Trial #8 RM/LL - 5mm	1	Z
5546-T-801	Tibial Augment Trial #8 LM/RL - 10mm	1	AA
5546-T-802	Tibial Augment Trial #8 RM/LL - 10mm	1	AB




Catalog #	Description	Qty	Ref. Letter
6543-7-522	IM Reamer - 22mm	1	A
6543-7-523	IM Reamer - 23mm	1	B
6543-7-524	IM Reamer - 24mm	1	C
6543-7-525	IM Reamer - 25mm	1	D
5566-T-022	Triathlon Stem Trial, 22 x 150mm	1	H
5566-T-023	Triathlon Stem Trial, 23 x 150mm	1	G
5566-T-024	Triathlon Stem Trial, 24 x 150mm	1	F
5566-T-025	Triathlon Stem Trial, 25 x 150mm	1	E
5565-T-022	Triathlon Stem Trial, 22 x 100mm	1	L
5565-T-023	Triathlon Stem Trial, 22 x 100mm	1	K
5565-T-024	Triathlon Stem Trial, 22 x 100mm	1	J
5565-T-025	Triathlon Stem Trial, 22 x 100mm	1	I




Catalog #	Description	Qty	Ref. Letter
5537-T-809	#8 9mm TS+ Insert Trial	1	A
5537-T-811	#8 11mm TS+ Insert Trial	1	B
5537-T-813	#8 13mm TS+ Insert Trial	1	C
5537-T-816	#8 16mm TS+ Insert Trial	1	D
5537-T-819	#8 19mm TS+ Insert Trial	1	E
5537-T-709	#7 9mm TS+ Insert Trial	1	F
5537-T-711	#7 11mm TS+ Insert Trial	1	G
5537-T-713	#7 13mm TS+ Insert Trial	1	H
5537-T-716	#7 16mm TS+ Insert Trial	1	I
5537-T-719	#7 19mm TS+ Insert Trial	1	J
5537-T-209	#2 9mm TS+ Insert Trial	1	K
5537-T-211	#2 11mm TS+ Insert Trial	1	L
5537-T-213	#2 13mm TS+ Insert Trial	1	M
5537-T-216	#2 16mm TS+ Insert Trial	1	N

Catalog #	Description	Qty	Ref. Letter
5537-T-219	#2 19mm TS+ Insert Trial	1	O
5537-T-222	#2 22mm TS+ Insert Trial	1	P
5537-T-225	#2 25mm TS+ Insert Trial	1	Q
5537-T-228	#2 28mm TS+ Insert Trial	1	R
5537-T-231	#2 31mm TS+ Insert Trial	1	S
5537-T-109	#1 9mm TS+ Insert Trial	1	T
5537-T-111	#1 11mm TS+ Insert Trial	1	U
5537-T-113	#1 13mm TS+ Insert Trial	1	V
5537-T-116	#1 16mm TS+ Insert Trial	1	W
5537-T-119	#1 19mm TS+ Insert Trial	1	X
5537-T-122	#1 22mm TS+ Insert Trial	1	Y
5537-T-225	#1 25mm TS+ Insert Trial	1	Z
5537-T-128	#1 28mm TS+ Insert Trial	1	AA
5537-T-131	#1 31mm TS+ Insert Trial	1	AB




Catalog #	Description	Qty	Ref. Letter
5537-T-722	#7 22mm TS+ Insert Trial	1	A
5537-T-725	#7 25mm TS+ Insert Trial	1	B
5537-T-728	#7 28mm TS+ Insert Trial	1	C
5537-T-731	#7 31mm TS+ Insert Trial	1	D
5537-T-822	#8 22mm TS+ Insert Trial	1	E
5537-T-825	#8 25mm TS+ Insert Trial	1	F
5537-T-828	#8 28mm TS+ Insert Trial	1	G
5537-T-831	#8 31mm TS+ Insert Trial	1	H




Catalog #	Description	Qty	Ref. Letter
6543-1-733	#3 Trial Cutting Guide	1	A
6543-1-734	#4 Trial Cutting Guide	1	B
6543-1-735	#5 Trial Cutting Guide	1	C
6543-1-736	#6 Trial Cutting Guide	1	D
5538-T-309	#3 TCG Insert 9mm	1	E
5538-T-311	#3 TCG Insert 11mm	1	F
5538-T-313	#3 TCG Insert 13mm	1	G
5538-T-409	#4 TCG Insert 9mm	1	H
5538-T-411	#4 TCG Insert 11mm	1	I
5538-T-413	#4 TCG Insert 13mm	1	J
5538-T-509	#5 TCG Insert 9mm	1	K
5538-T-511	#5 TCG Insert 11mm	1	L

Catalog #	Description	Qty	Ref. Letter
5538-T-513	#5 TCG Insert 13mm	1	M
5538-T-609	#6 TCG Insert 9mm	1	N
5538-T-611	#6 TCG Insert 11mm	1	O
5538-T-613	#6 TCG Insert 13mm	1	P
6543-1-740	Neutral TCG Valgus Adapter	1	Q
6543-1-742	2mm Offset TCG Valgus Adapter	1	R
6543-1-744	4mm Offset TCG Valgus Adapter	1	S
6543-1-746	6mm Offset TCG Valgus Adapter	1	T
6543-1-748	8mm Offset TCG Valgus Adapter	1	U
5100-3600	TCG Torque Handle	1	V
6543-4-820	TCG T-20 Torx Driver	1	W
6543-4-524	TCG Reamer - 24mm	1	X




Catalog #	Description	Qty	Ref. Letter
5538-T-316	#3 TCG Insert 16mm	1	A
5538-T-319	#3 TCG Insert 19mm	1	B
5538-T-322	#3 TCG Insert 22mm	1	C
5538-T-325	#3 TCG Insert 25mm	1	D
5538-T-416	#4 TCG Insert 16mm	1	E
5538-T-419	#4 TCG Insert 19mm	1	F
5538-T-422	#4 TCG Insert 22mm	1	G
5538-T-425	#4 TCG Insert 25mm	1	H
5538-T-516	#5 TCG Insert 16mm	1	I
5538-T-519	#5 TCG Insert 19mm	1	J
5538-T-522	#5 TCG Insert 22mm	1	K
5538-T-525	#5 TCG Insert 25mm	1	L
5538-T-616	#6 TCG Insert 16mm	1	M
5538-T-619	#6 TCG Insert 19mm	1	N
5538-T-622	#6 TCG Insert 22mm	1	O
5538-T-625	#6 TCG Insert 25mm	1	P




Catalog #	Description	Qty	Ref. Letter
6543-1-731	#1 Trial Cutting Guide	1	A
6543-1-732	#2 Trial Cutting Guide	1	B
6543-1-737	#7 Trial Cutting Guide	1	C
6543-1-738	#8 Trial Cutting Guide	1	D
6543-1-740	Neutral TCG Valgus Adapter	1	E
6543-1-742	2mm Offset TCG Valgus Adapter	1	F
6543-1-744	4mm Offset TCG Valgus Adapter	1	G
6543-1-746	6mm Offset TCG Valgus Adapter	1	H
6543-1-748	8mm Offset TCG Valgus Adapter	1	I
5100-3600	TCG Torque Handle	1	J
6543-4-820	TCG T-20 Torx Driver	1	K
6543-4-524	TCG Reamer - 24mm	1	L

Catalog #	Description	Qty	Ref. Letter
5538-T-109	#1 TCG Insert 9mm	1	M
5538-T-111	#1 TCG Insert 11mm	1	N
5538-T-113	#1 TCG Insert 13mm	1	O
5538-T-209	#2 TCG Insert 9mm	1	P
5538-T-211	#2 TCG Insert 11mm	1	Q
5538-T-213	#2 TCG Insert 13mm	1	R
5538-T-709	#7 TCG Insert 9mm	1	S
5538-T-711	#7 TCG Insert 11mm	1	T
5538-T-713	#7 TCG Insert 13mm	1	U
5538-T-809	#8 TCG Insert 9mm	1	V
5538-T-811	#8 TCG Insert 11mm	1	W
5538-T-813	#8 TCG Insert 13mm	1	X



Catalog #	Description	Qty	Ref. Letter
5538-T-116	#1 TCG Insert 16mm	1	A
5538-T-119	#1 TCG Insert 19mm	1	B
5538-T-122	#1 TCG Insert 22mm	1	C
5538-T-125	#1 TCG Insert 25mm	1	D
5538-T-216	#2 TCG Insert 16mm	1	E
5538-T-219	#2 TCG Insert 19mm	1	F
5538-T-222	#2 TCG Insert 22mm	1	G
5538-T-225	#2 TCG Insert 25mm	1	H
5538-T-716	#7 TCG Insert 16mm	1	I
5538-T-719	#7 TCG Insert 19mm	1	J
5538-T-722	#7 TCG Insert 22mm	1	K
5538-T-816	#8 TCG Insert 16mm	1	L
5538-T-819	#8 TCG Insert 19mm	1	M
5538-T-822	#8 TCG Insert 22mm	1	N



Catalog #	Description	Qty	Ref. Letter
5538-T-128	#1 TCG Insert 28mm	1	A
5538-T-131	#1 TCG Insert 31mm	1	B
5538-T-228	#2 TCG Insert 28mm	1	C
5538-T-231	#2 TCG Insert 31mm	1	D
5538-T-328	#3 TCG Insert 28mm	1	E
5538-T-331	#3 TCG Insert 31mm	1	F
5538-T-428	#4 TCG Insert 28mm	1	G
5538-T-431	#4 TCG Insert 31mm	1	H
5538-T-528	#5 TCG Insert 28mm	1	I

Catalog #	Description	Qty	Ref. Letter
5538-T-531	#5 TCG Insert 31mm	1	J
5538-T-628	#6 TCG Insert 28mm	1	K
5538-T-631	#6 TCG Insert 31mm	1	L
5538-T-725	#7 TCG Insert 25mm	1	M
5538-T-728	#7 TCG Insert 28mm	1	N
5538-T-731	#7 TCG Insert 31mm	1	O
5538-T-825	#8 TCG Insert 25mm	1	P
5538-T-828	#8 TCG Insert 28mm	1	Q
5538-T-831	#8 TCG Insert 31mm	1	R

A surgeon must always rely on his or her own professional clinical judgment when deciding to use which products and/or techniques on individual patients. Stryker is not dispensing medical advice and recommends that surgeons be trained in orthopaedic surgeries before performing any surgeries.

The information presented is intended to demonstrate the breadth of Stryker product offerings. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker and Triathlon. All other trademarks are trademarks of their respective owners or holders.

Literature Number: **LTTSPC**
MS/GS 5c 12/07

Copyright © 2007 Stryker
Printed in USA



Orthopaedics

325 Corporate Drive
Mahwah, NJ 07430
t: 201.831.5000

www.stryker.com